

Product datasheet for **SC200866**

PRICKLE3 (NM_006150) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: PRICKLE3

Synonyms: LMO6; LOAM; Pk3

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PSI00062)

ACCN: NM_006150

Insert Size: 901 bp

Insert Sequence: >SC200866 3'UTR clone of NM_006150
The sequence shown below is from the reference sequence of NM_006150. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GCCCGAGACAAGAACTGCATCGTGGCTGAAGGCAGGCCGTCTGGAGGGGGCTCCATTCTCCAGTCAG
AGTAGATGATGAGGCCCATGCCCTCACCCACGCCCGCCCTACAACCTAAGTCATAAATCCTCTT
CCTCCCTCCTTTATCAGTTTGTCTATTTTGTCTGTTTGTAGAGACTTGGTCTACCTATTGCCCA
GGCTGGAGTGCAGTGGTGTGATCAGGGCTCACTGCAGCCTCAAAAGTTTGGGCTCAAGCTATCCTCCC
ACCTCAGCCTCCCAAGTAGCTGGGACTCCAGGCTCGTGTGAGCGTTCCAGCTAATTTTAAATTTT
GGTAGAGATGGGGTTTGCCATGTTGCCAGGCTGGTCTAGAACTCCTGGGCTCAAGCAATCCTCCCGC
AGCCTCCCAAGGACTGGGATTACAGGCATGAGCCACCGCACCAGGCTCTCGTTTGTATTATTGAAGT
AATTTAATATTGTGTAAACAGGCTGGCTTCTTTGCCTCCTTGCCCAAGCTGTTTCCACTCCATTT
CCCCCTCCTGATGCTTGAAATTTACCCACATTTGATGAAAGGATTACTTTATTCAATTATCGTATATTTT
AACAATGTTTATTATTCATTTATCCCTCTATAGAACCACACCCACAGGAGATTATTTGGAGTGG
GTCCCAACCTAGGGCTGGACTCTGAAATCTAACTCCCACTTCCCTCATTTTGTGACTTAGGTGGGG
CATGGTTTCAGTCAGAACTGGTGTCTCTATTGGATCGTGCAGAAGGAGGACCTAGGCACACACATATGG
TGGCCACACCCAGGAGGTTGATTGGCAGGCTGGAAGACAAAAGTCTCCAATAAAGGCACTTTTACCT
CAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_006150.5</u>
Summary:	LIM domain only 6 is a three LIM domain-containing protein. The LIM domain is a cysteine-rich sequence motif that binds zinc atoms to form a specific protein-binding interface for protein-protein interactions. [provided by RefSeq, Jul 2008]
Locus ID:	4007
MW:	33.4