

Product datasheet for **SC200900**

ZNF342 (ZNF296) (NM_145288) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF342
Synonyms:	ZFP296; ZNF342
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_145288
Insert Size:	131 bp
Insert Sequence:	>SC200900 3'UTR clone of NM_145288 The sequence shown below is from the reference sequence of NM_145288. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC AAGCACCCCTGAGCGCGCCGGCGAGGCCGAGGCCAGGCCAGGAAAGCCCCCTCACTGTCCCTGGTACCGCTGC CAACACCCATTGACCTCTCGTTTTTGCCCGCCTTCTCCAAGTAAATTTCCCTTTTATTTA ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.



RefSeq: NM_145288.3

Summary: May be a transcriptional corepressor with KLF4.[UniProtKB/Swiss-Prot Function]

Locus ID: 162979

MW: 4.9