

Product datasheet for **SC200936**

ASNSD1 (NM_019048) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ASNSD1
Synonyms:	NBLA00058; NS3TP1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_019048
Insert Size:	135 bp
Insert Sequence:	>SC200936 3'UTR clone of NM_019048 The sequence shown below is from the reference sequence of NM_019048. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CTTTCTATTGAAAAGGAGACTAAATTGAAATGTGATTACAAATGTAACAATATAAAAAATAAGTTTTTAT ATAATTATATAAAAGTAAGATACTCTGCTGCTTTACTATTGTATAACATAGTAGTTTTAAAGTTCA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.



RefSeq: NM_019048.4

Locus ID: 54529

MW: 5.3