

Product datasheet for **SC200948**

FKBP2 (NM_004470) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FKBP2 (NM_004470) Human 3' UTR Clone
Symbol:	FKBP2
Synonyms:	FKBP-13; FKBP13; PPIase
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004470
Insert Size:	133 bp
Insert Sequence:	>SC200948 3'UTR clone of NM_004470 The sequence shown below is from the reference sequence of NM_004470. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC CTCAAAATAGAGCGACGAACTGAGCTGT AA CCAGACTGGGGAGGGGCAGGGGGAGAGGCCCCCATCAGG GACCAGACTGTTCAAAAAAAAAACAAAAACAAAAACAAAAACAAAAACACTTAAAGCCCAA ACGCGT AAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_004470.4</u>



[View online »](#)

Summary:

The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin. It is thought to function as an ER chaperone and may also act as a component of membrane cytoskeletal scaffolds. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Sep 2008]

Locus ID:

2286

MW:

4.9