

## Product datasheet for **SC200968**

### CUTA (NM\_015921) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	CUTA (NM_015921) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CUTA
Synonyms:	ACHAP; C6orf82
ACCN:	NM_015921
Insert Size:	238 bp
Insert Sequence:	>SC200968 3'UTR clone of NM_015921 The sequence shown below is from the reference sequence of NM_015921. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCA <b>ACGATCGCC</b> GTTTCTGACTCTATCACAGTCTGCCA <b>TGAT</b> GAGCCCTGTTCTGCTCATCATGAAGATCCCCGCGATA CTTCAACGCCTTCTGACTTCCAGGTGATGACTGGGCCCAATAAATCCCGTCTTTGGGTCTCTGCC TTTTGCCTGTTTTTGTGATGTATGCAAATTTGTTGGGAGAATATAAGTGCTTCTAGTCCCTGTATTCT CTACATCTTTCAGAATGACTCAGAGCTGCTA <b>ACGCGT</b> AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<a href="#">NM_015921.3</a>
Summary:	May form part of a complex of membrane proteins attached to acetylcholinesterase (AChE). [UniProtKB/Swiss-Prot Function]
Locus ID:	51596



[View online »](#)

MW: 8.8