

Product datasheet for **SC200971**

CUTA (NM_001014840) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CUTA (NM_001014840) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CUTA
Synonyms:	ACHAP; C6orf82
ACCN:	NM_001014840
Insert Size:	238 bp
Insert Sequence:	>SC200971 3'UTR clone of NM_001014840 The sequence shown below is from the reference sequence of NM_001014840. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCA ACGATCGCC GTTTCTGACTCTATCACAGTCCTGCCAT GT AGAGCCCTGTTCTGCTCATCATGAAGATCCCCGCGATA CTTCAACGCCCTTCTGACTTCCAGGTGATGACTGGGCCCAATAAATCCCGTCTTTGGGTCTCTCTGCC TTTTGCCTGTTTTTGTATGTATGCAAATTTGTTGGGAGAATATAAGTGCTTCTAGTCCCTGTATTCT CTACATCTTTCAGAATGACTCAGAGCTGCTA ACGCGT AAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001014840.2</u>
Summary:	May form part of a complex of membrane proteins attached to acetylcholinesterase (AChE). [UniProtKB/Swiss-Prot Function]
Locus ID:	51596



[View online »](#)

MW: 8.8