

Product datasheet for **SC201035**

FAM71B (NM_130899) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: FAM71B

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PS100062)

ACCN: NM_130899

Insert Size: 631 bp

Insert Sequence: >SC201035 3'UTR clone of NM_130899
The sequence shown below is from the reference sequence of NM_130899. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATGATCGTCTTTGAAACCAATCCATTAAAGAGGAGCCAGGGCTACAACGCCCCAGGGATCTGGAAT
ATCCCTAGAGGAGCCCATCCCAGCTTTCTTTACTGCAGTTCAAATAATAAGGAAATCGTACAACCTCCA
AAAGAAATACAATCTCCTGCTTGAATTTTTCTGCCCATAGCCTGGTGGTGACCTCACAGTGGCTTCTGT
GATGGTGGCTATGCTACAGCTGGAGACCCAGCCTCCAGGCAAGTAGAGCTCCACCACCCGCCACGGC
CAGCAGAGCCAACGTGAGTGTTTCAGCTTTGCCTCAATCCCTTTTGCTTCCCGAGACTGTCAAAGCCAC
AGACTAAGAGCTGGGATGCATGGCCACCACCATGTCCCTCGCCCATTTCTCCTCTTTTCCCTGAATGCC
CCTGAACATGGACATTCAGTCCAGTTCTTAGCCTATGCAGAGACAGGATTCACATCCCTCAGGGGATA
TTCAGGGTAGGAAGGAGCAATACAGCAGTATTTGTAGGCAATTAGAATAGGAAATCTATTCTTCTAGAT
CAGGAATTGGCAAATCATGTTCTGAGGTCCAAATCTTACCAGCTGCCTGGATTTGTAATAAAGTTTAA
TTGGAACATA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_130899.3</u>
Summary:	May be involved in RNA biogenesis.[UniProtKB/Swiss-Prot Function]
Locus ID:	153745
MW:	22.9