

## Product datasheet for SC201090

### UBL5 (NM\_001048241) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	UBL5 (NM_001048241) Human 3' UTR Clone
Symbol:	UBL5
Synonyms:	HUB1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001048241
Insert Size:	137 bp
Insert Sequence:	<p>&gt;SC201090 3'UTR clone of NM_001048241</p> <p>The sequence shown below is from the reference sequence of NM_001048241. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GGGATGAACCTGGAGCTTTATTATCAATAGATGAGAATCCTCATCTTCTGCCCCGCTTTCCTCTCCCA TCCTCATCCCCACACTGGGATAGATGCTTGTGTAATAAACTCACCTTAATAAAGACTTAGATGTTG ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG           </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<a href="#">NM_001048241.3</a>


[View online »](#)

**Summary:**

This gene encodes a member of a group of proteins similar to ubiquitin. The encoded protein is not thought to degrade proteins like ubiquitin but to affect their function through being bound to target proteins by an isopeptide bond. The gene product has been studied as a link to predisposition to obesity based on its expression in Psammomys obesus, the fat sand rat, which is an animal model for obesity studies. Variation in this gene was found to be significantly associated with some metabolic traits (PMID: 15331561) but not associated with childhood obesity (PMID: 19189687). Pseudogenes of this gene are located on chromosomes 3, 5 and 17. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jan 2013]

**Locus ID:**

59286

**MW:**

5.2