

Product datasheet for **SC201107**

LRRC36 (NM_001161575) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: LRRC36

Synonyms: RORBP70; XLHSRF2

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_001161575

Insert Size: 145 bp

Insert Sequence: >SC201107 3'UTR clone of NM_001161575
The sequence shown below is from the reference sequence of NM_001161575. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTGGATGCTGCCCCAGAGCCTGGCTTAAGAGCTAGCATGGAACACACACAGCTTCCCTGGTCCAC
AGAGGCTCTACCGCCATTGCCACCATGTTGGTGTGTACTCACAAGATTAAGAAAGAAATGTATT
CTGATTA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001161575.2

Locus ID: 55282

MW: 5.4