

Product datasheet for **SC201132**

LAIR2 (NM_002288) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LAIR2 (NM_002288) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	LAIR2
Synonyms:	CD306
ACCN:	NM_002288
Insert Size:	148 bp
Insert Sequence:	>SC201132 3'UTR clone of NM_002288 The sequence shown below is from the reference sequence of NM_002288. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC ACTGAAGCCTCCGGATTTGATGCACCATGAAATGAGGAGAAATGGCCTCCCGTCTTGTGAACTTCAATGG GGAGAAATAATTAGAATGAGCAATAGAAATGCACAGATGCCTATACATACATATACAAATAAAAAGATA CGATTGCAA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_002288.6</u>



[View online »](#)

Summary: The protein encoded by this gene is a member of the immunoglobulin superfamily. It was identified by its similarity to leukocyte-associated immunoglobulin-like receptor 1, a membrane-bound receptor that modulates innate immune response. The protein encoded by this locus is a soluble receptor that may play roles in both inhibition of collagen-induced platelet aggregation and vessel formation during placental implantation. This gene maps to a region of 19q13.4, termed the leukocyte receptor cluster, which contains 29 genes in the immunoglobulin superfamily. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Sep 2013]

Locus ID: 3904

MW: 5.7