

Product datasheet for **SC201164**

ALGIL2 (NM_001136152) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ALGIL2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001136152
Insert Size:	165 bp
Insert Sequence:	>SC201164 3'UTR clone of NM_001136152 The sequence shown below is from the reference sequence of NM_001136152. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TTTGAGATGCTTTTCTCAAATTTCTGATCCTGAAGCAAGCTAAACCAGTTCCGGAAGAACCTGGC GGAGTCGCAGCAGCTCTGATGGGATGAGAGCTGGGTGCAGACTGTGCTCCCTTTGGTTATGGACACATA ACTCCTGGGCCAGAGGCTAAAACCCCA ACGCGTAAGCGGCCCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.



RefSeq: NM_001136152.1

Summary: Putative glycosyltransferase.[UniProtKB/Swiss-Prot Function]

Locus ID: 644974

MW: 5.6