

Product datasheet for **SC201171**

MZT2B (NM_025029) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	MZT2B
Synonyms:	FAM128B; MOZART2B
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_025029
Insert Size:	131 bp
Insert Sequence:	>SC201171 3'UTR clone of NM_025029 The sequence shown below is from the reference sequence of NM_025029. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GGGAAGAGCCCTACACGGGGCAGCACCTAGGATGGGGCAGAGACTTGTTCATCTTTGTCCCAGCAAA GGCTACATGTTACCTCCTTCAATTGATAATAAACCTTTCTGAGATGCAGAGGGTCCAGGTCA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.



RefSeq: NM_025029.5

Locus ID: 80097

MW: 4.5