

Product datasheet for **SC201187**

LRTM1 (NM_020678) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: LRTM1

Synonyms: HT017

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PSI00062)

ACCN: NM_020678

Insert Size: 258 bp

Insert Sequence: >SC201187 3'UTR clone of NM_020678
The sequence shown below is from the reference sequence of NM_020678. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AAAGAGCGATTTGACAGCTCACCAGCCAGAGCTTTTGTCTCAAATAGGATTGGTCATTGCAGGCCAG
AAGATAGTGTCTGAGTAGGGCTGATGTGTTCTGTAGTCTGATTTGCTTTTGCCAAAAGACAAAAA
AAAAAAAAAAAAAAAAAGAAAAAAAAAGATGTGTCTGTAAAAATATTTCTGGGGAGGAATTTATGT
TGCGATTTTAAATACTATAGGTAAATAAAGATGATACTTGTGCTCTGG
ACGCGTAAGCGGCCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_020678.4</u>
Locus ID:	57408
MW:	9.8