

Product datasheet for SC201188

THAP7 (NM_030573) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	THAP7 (NM_030573) Human 3' UTR Clone
Symbol:	THAP7
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_030573
Insert Size:	806 bp
Insert Sequence:	>SC201188 3'UTR clone of NM_030573 The sequence shown below is from the reference sequence of NM_030573. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GCCATGCAGCTGAGCAGCAGCATGGCCTGAGGGGCTGCTGGACTGACCGAGGGGCTGCCAGCAAGACT
GCAGCCTCTTCTCCCTCAGATCCCACCAGACCCACAGGTGCCATAATAAGCGGATTCTAGACGGAG
GGCCTGCTGTCTGGCTCAGCTGGGGCTGGAGCCCCTCTGAGTACTAAGCTTGCTGCCCCAGACAGTGG
CTGCAGTCTCTAGAGGGCTCAGTCCAAGCTCTTTGCCAAGGAGAACAGGGAGCCTGGCTCTCACCCTC
AGAGCTACTGGTAGACAGGTGGCAGTCAGCCCTGCTAGCGCAGCTCCGTGCCCCACAGGGTTCTCTGAA
GGGGCAGGATGGTGTGAGTAGGGGCCAGGAAGACAGCTCCAGGCGCCCTTTCATAGGATCTGAAACAG
CCCTGACTTAAACCAGCCAGCCCTTTGAGAACACTTCTAAAACCCAGTCTAGATCCGGTTCTCAAGCA
GGAGCATGACAGGCTTCTGTCCACAGGCCCCATCAGCCTGGTGGGTAAGCTGGACAACCCTCTCCCCGC
CCCCTTCACTACCATCTCTATTTGAAAAACGGGGATGATGATGACAATCATCAACCTCCTGGGTCTGT
CAGCAGGTAGTAAAGCCCCAACAAGCTAGTCCCTGGAAATCTAGGAAGAGGTGGCCAGAATGGAGAC
AAGCACTCAGGTCTGGCTCAGCATCTTTTCTTCAACTCCCAGGAACCTGGAAGGCCTTGCTGCAG
CTGCAGTTCTCAAACAGTCCGTCTTGGGGTGACTCACTGCAAGTAT
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_030573.3
Summary:	Chromatin-associated, histone tail-binding protein that represses transcription via recruitment of HDAC3 and nuclear hormone receptor corepressors.[UniProtKB/Swiss-Prot Function]
Locus ID:	80764
MW:	28.4