

Product datasheet for **SC201214**

RBM6 (NM_005777) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	RBM6
Synonyms:	3G2; DEF-3; DEF3; gl6; HLC-11; NY-LU-12
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_005777
Insert Size:	146 bp
Insert Sequence:	>SC201214 3'UTR clone of NM_005777 The sequence shown below is from the reference sequence of NM_005777. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC ATGTTTGCTCGATATAAAGAACTCGATTAAAGAAAGGAGACAAGTTCCATGGGATACAACCTCCCTCTTG TTTTGTTTGCTCTCCTTTTCTTTTGTACTGTTCTTGCTGCTAGAACTTTTTTAAATAAACTTTTTTT CAATGTGA ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_005777.3</u>
Summary:	Specifically binds poly(G) RNA homopolymers in vitro.[UniProtKB/Swiss-Prot Function]
Locus ID:	10180
MW:	5.9