

Product datasheet for **SC201266**

ARSA (NM_001085427) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ARSA (NM_001085427) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ARSA
Synonyms:	ASA; MLD
ACCN:	NM_001085427
Insert Size:	2000 bp



[View online »](#)

Summary: The protein encoded by this gene hydrolyzes cerebroside sulfate to cerebroside and sulfate. Defects in this gene lead to metachromatic leucodystrophy (MLD), a progressive demyelination disease which results in a variety of neurological symptoms and ultimately death. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Dec 2010]

Locus ID: 410

MW: 72.5