

Product datasheet for **SC201279**

Aurora C (AURKC) (NM_001015878) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Aurora C (AURKC) (NM_001015878) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	AURKC
Synonyms:	AIE2; AIK3; ARK3; AurC; aurora-C; HEL-S-90; SPGF5; STK13
ACCN:	NM_001015878
Insert Size:	161 bp
Insert Sequence:	>SC201279 3'UTR clone of NM_001015878 The sequence shown below is from the reference sequence of NM_001015878. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CTGCCTCCCTGTGCTCAGATGGCTTCTGAGCCCTGTCTGCCTCTGTTCCCTTTGTGTGTTCAGGGA GCTCTCCTGGCTCTGCCACCTCATTGTCTTTATTTTTTCTCTTTAAGATGTAAGATGCTAATTAAT AAAAGCTGAATCATTTCATACCA ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001015878.2</u>



[View online »](#)

Summary: This gene encodes a member of the Aurora subfamily of serine/threonine protein kinases. The encoded protein is a chromosomal passenger protein that forms complexes with Aurora-B and inner centromere proteins and may play a role in organizing microtubules in relation to centrosome/spindle function during mitosis. This gene is overexpressed in several cancer cell lines, suggesting an involvement in oncogenic signal transduction. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008]

Locus ID: 6795

MW: 5.7