

Product datasheet for **SC201291**

ALAS2 (NM_001037968) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ALAS2 (NM_001037968) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ALAS2
Synonyms:	ALAS-E; ALASE; ANH1; ASB; SIDBA1; XLDPP; XLEPP; XLSA
ACCN:	NM_001037968
Insert Size:	155 bp
Insert Sequence:	>SC201291 3'UTR clone of NM_001037968 The sequence shown below is from the reference sequence of NM_001037968. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GGGCCCCAGTATGTCACCACCTATGCC TGAAGCCAGCTGCCTAGGATTCACACCCACCTGCGCTTC ACTTGGGTCCAGGCCTACTCCTGTCTTCTGCTTTGTTGTGTGCCTCTAGCTGAATTGAGCCTAAAAATA AAGCACAAACCACAGCA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001037968.4</u>



[View online »](#)

Summary: The product of this gene specifies an erythroid-specific mitochondrially located enzyme. The encoded protein catalyzes the first step in the heme biosynthetic pathway. Defects in this gene cause X-linked pyridoxine-responsive sideroblastic anemia. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Locus ID: 212

MW: 5.7