

Product datasheet for **SC201315**

DHPS (NM_013407) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DHPS (NM_013407) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	DHPS
Synonyms:	DHS; MIG13
ACCN:	NM_013407
Insert Size:	174 bp
Insert Sequence:	>SC201315 3'UTR clone of NM_013407 The sequence shown below is from the reference sequence of NM_013407. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GCCTTCATGCATGAGAAGAACGAGGACTGAGCGGCTGCGGTCCCAGGAAGGTCTTACCCCTCTTCTAT TTATTAATTTGCAGACCCAGCCCTCCCTACTTTTTGGTCAGCTACGTCTCTAGAATAAGATGGTATC TGAAGTCCTTAAAAAAAAAAAAAAAAAAAAAAAAAAAA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_013407.1</u>



[View online »](#)

Summary: This gene encodes a protein that is required for the formation of hypusine, a unique amino acid formed by the posttranslational modification of only one protein, eukaryotic translation initiation factor 5A. The encoded protein catalyzes the first step in hypusine formation by transferring the butylamine moiety of spermidine to a specific lysine residue of the eukaryotic translation initiation factor 5A precursor, forming an intermediate deoxyhypusine residue. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, May 2011]

Locus ID: 1725

MW: 6.6