

Product datasheet for SC201338

DEPDC7 (NM_001077242) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: DEPDC7

Synonyms: dJ85M6.4; TR2

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PSI00062)

ACCN: NM_001077242

Insert Size: 157 bp

Insert Sequence: >SC201338 3'UTR clone of NM_001077242
The sequence shown below is from the reference sequence of NM_001077242. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GACATCTTTATTGAGCATTTTGGAGACTGAGTTTTTAATATCTGTATATAAGTTGTGTATTTAAGAAT
AAATTATGTATCCTAAATATCCAATCACATTTGTAAGCGTGGAAGCTCTAAATTTGAAACTGTACTTAA
TAAAAATTTTTTTGTATAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001077242.2](#)

Locus ID: 91614

MW: 6.2