

Product datasheet for **SC201349**

MAP3K19 (NM_025052) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: MAP3K19 (NM_025052) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: MAP3K19
Synonyms: RCK; YSK4
ACCN: NM_025052
Insert Size: 389 bp
Insert Sequence: >SC201349 3'UTR clone of NM_025052

The sequence shown below is from the reference sequence of NM_025052. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGCACTCCTTCTTGGAGAGAAGTCACATGATATACATCAAGACTTCTTCCCAGTTCCTACTGCAGATG
CTCCCTTGCTTAATTGTGGGAATGATGGCTAAGGGATCTTTGTTTCCCCTGAAAATTCAGTCTAAC
CCAGTTTAAGCAGATCCTATGGAGTCATTAAGTAAAGTTGCAGTTACATATTAGCTCCTCAAGTGTC
AGACATTATTACTCATAGTATCAGAAAACATGTTCTTAATAACAACAAAAAACTATTTTCAGTGTTTACA
GTTTTGATTGTCCAGGAACACTACATTCTCTATTGTTTTATATGACATTTCTTTTTATTTTTGGCCTGTCC
TGTCATTTTTAATGTTGTTAGTTTAAAATAAATTGTAACAAACAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_025052.5](#)

Locus ID: 80122



[View online >](#)

MW: 14.8