

Product datasheet for **SC201350**

MAP3K19 (NM_001018046) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MAP3K19 (NM_001018046) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	MAP3K19
Synonyms:	RCK; YSK4
ACCN:	NM_001018046
Insert Size:	389 bp
Insert Sequence:	>SC201350 3'UTR clone of NM_001018046 The sequence shown below is from the reference sequence of NM_001018046. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC AAGCACTCCTTCTGGAGAGAAGTCACT GA ATATACATCAAGACTTCTCCAGTTCCTACTGCAGATG CTCCCTTGCTTAATTGTGGGAATGATGGCTAAGGGATCTTTGTTCCCACTGAAAATTCAGTCTAAC CCAGTTTAAGCAGATCCTATGGAGTCATTAAGTAAAGTTGCAGTTACATATTAGCTCCTCAAGTGTC AGACATTACTCATAGTATCAGAAAACATGTTCTTAATAACAACAAAAAACTATTTTCAGTGTTTACA GTTTTGATTGTCCAGGAACACTATTCTCTATTGTTTTATATGACATTTCTTTTTATTTTGGCCTGTCC TGTCATTTTTAATGTTGTTAGTTTAAAATAAATTGTA AAAACAA ACGCGT AAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001018046.3</u>
Locus ID:	80122



MW: 14.8