

Product datasheet for **SC201354**

ERGIC3 (NM_198398) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ERGIC3
Synonyms:	C2orf47; C20orf47; CGI-54; dJ477O4.2; Erv46; NY-BR-84; PRO0989; SDBCAG84
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_198398
Insert Size:	160 bp
Insert Sequence:	>SC201354 3'UTR clone of NM_198398 The sequence shown below is from the reference sequence of NM_198398. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC AAGAAAATTGATCTAGGGAAGACAACGAGTACACCTCGGTGCTTCTGTCTCTTCTCCCTGG CCTGTGGTTGTCCCCAGCCTCTGCCACCCTCCACCTCCTCGGTGACGCCAGCCAGGTTGATAAT CTATTGATTGATTGTAGTAGTA ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACGCCGCTTCTATGAAAG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_198398.2</u>
Summary:	Possible role in transport between endoplasmic reticulum and Golgi.[UniProtKB/Swiss-Prot Function]
Locus ID:	51614
MW:	5.6