

Product datasheet for **SC201362**

EVX1 (NM_001989) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	EVX1 (NM_001989) Human 3' UTR Clone
Symbol:	EVX1
Synonyms:	EVX-1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001989
Insert Size:	1435 bp



[View online »](#)

Insert Sequence: >SC201362 3'UTR clone of NM_001989
The sequence shown below is from the reference sequence of NM_001989. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGAGGGAGGAGGTGCCCTCACTAGATAAGGGGCCGCGGCTGGCTGCCGCTCCATGACGCCCGTGG
GGTCACCCCCCGCCCCGGGACTCAGCCAGCCTCGCTCCTCGCTCCTCGCTCCTCGCCCCCTAGGACGCC
AAGGGGAAAGGAGAGGGCGGAAAGGACCAGCGGGATCCGGCCGCAAGAATTGGAAGCCTAGGAAGT
GGCGGTGGCTGGCGCTTTGGGAGCAGGAGTGGGGATAGGGAAGCAGAGCTTGAGAGACCTTCTCCG
GGGCAGCTCCGACCCACCGCCCCCACCAGGTCGAGGCTGTAGCTCCAAAGCTAAACAAAATTAG
CAGCAACAGCAACCAATATCCAGTCCCTCGGCCCTCGGCCCTCACCTCCACCTCACACTCCCTTCT
CACCGGGCCCCCTCTCCCAGCCAAGGCCAAGCACTGGAAGGGAAATTGCTGTCTCTGTAACAAAA
TGCTGTGTATGCAGAGCAGGTAGAGATTAATCTTTGCCAGCTTTTCCAAGGCATGACAAGGGGCTGGTG
GATGGCAACATACCAGTCATTTGGAGGAGAGAGTGAGAGATGATTACTACAGGAGAAATCCAGCCCC
TTGGCATGGGACCTGGAGCCTCGACTACACAGCATCTTCTGGGTCTGGCCTCTGCCAGCACCTGATCTC
TTTCTCATTCCAGCTTTGTGACACTTCTCAACTTGGCGCTCCATCTCTCCCTGCCCCCACTTTTTTG
TTGGCCAGGGAGGCTGCAGATGCCCCAGGAGCCCTTTGCCGCTTCTATGAGGCCAAGCCTTTTTTCCT
GGGCCCAGCACACCCTGATTAGCAAGTGATGTGTGCGAGGAGGGTTTGTGAATGTTGAATGTGAAT
AATGATCACCATGGAGCTGGCCACTGACCCAGAGCTGAGCTGTTAACAAGGCGCCAGGGAAGAGCTT
AGGGAGTGGGAACTTCACCTCCCTCTCTCGGTATCTGGCGGTAATTAAGAGCAATTTTCATCCTTTGC
TTGTTTACCTTCACTTCACCAGGAATTTCTGGCCCTACCCTTTGCATTGGGTATTTTACAATTTCTC
TCATTTTCTTCCCAAGCTACCACTGGAGCTTGACTTTTCAGATACCAGTGGGAGCCTTCTGTCCCTTTTG
GGGACCTGTCTGTGGCCTCCACCAGGGTTTGTAGAGCCACTCCCAATCCTCACTCCCACACTCAT
CCTTGACGCCAGTTTTTGAGGAAGAGGAGAACGTGAACCCCAATGCAAGCTTCACCCTGACTGAGAGG
GAGTGGTTCTTCTGTAGGGAATGAATTTGGTTTGATTGGGGTTTTCTTTGAAGCCCAAGAAGCTTG
CTGTTATGATTGTTAACCATATTGCAATAAAGCTGGACATAATTCTCACTTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001989.5](#)

Summary: This gene encodes a member of the even-skipped homeobox family characterized by the presence of a homeodomain closely related to the Drosophila even-skipped (eve) segmentation gene of the pair-rule class. The encoded protein may play an important role as a transcriptional repressor during embryogenesis. [provided by RefSeq, Jul 2008]

Locus ID: 2128

MW: 53