

Product datasheet for SC201400

C11orf71 (NM_019021) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: C11orf71 (NM_019021) Human 3' UTR Clone
Symbol: C11orf71
Synonyms: URLC7
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_019021
Insert Size: 152 bp
Insert Sequence: >SC201400 3'UTR clone of NM_019021
 The sequence shown below is from the reference sequence of NM_019021. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TCAAATATTTTGAATTTTAGTAATAAAGTACTGATTTGTGATATAATTACATTTATATTATTTAAGAT
AAAAATGGATATGAAATTGTATGTATTATATGATCTTACCTATATGAATGAAATATGCATTAATAAAC
TAGAATTGCCAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq: [NM_019021.4](#)
Locus ID: 54494


[View online »](#)

MW:

6.1