

Product datasheet for **SC201456**

Cytokeratin 17 (KRT17) (NM_000422) Human 3' UTR Clone

Product data:

| | |
|--------------------|--|
| Product Type: | 3' UTR Clones |
| Product Name: | Cytokeratin 17 (KRT17) (NM_000422) Human 3' UTR Clone |
| Vector: | pMirTarget (PS100062) |
| Symbol: | KRT17 |
| Synonyms: | 39.1; CK-17; K17; PC; PC2; PCHC1 |
| ACCN: | NM_000422 |
| Insert Size: | 182 bp |
| Insert Sequence: | >SC201456 3'UTR clone of NM_000422 The sequence shown below is from the reference sequence of NM_000422. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CGCGAGCAGGTCCACCAGACCACCGCTGAGGACTCAGCTACCCCGCCGCGCCACCCAGGAGGCAGGGA GGCAGCCGCCCATCTGCCACAGTCTCCGGCCTCTCCAGCCTCAGCCCCTGCTTCAGTCCCTTCCC CATGCTTCCCTTGCTGATGACAATAAAGCTTGTTGACTCAGCTA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG |
| Restriction Sites: | Sgfl-MluI |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | NM_000422.3 |
| Summary: | This gene encodes the type I intermediate filament chain keratin 17, expressed in nail bed, hair follicle, sebaceous glands, and other epidermal appendages. Mutations in this gene lead to Jackson-Lawler type pachyonychia congenita and steatocystoma multiplex. [provided by RefSeq, Aug 2008] |



[View online »](#)

Locus ID: 3872

MW: 6.5