

Product datasheet for **SC201477**

DAB2IP (NM_138709) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: DAB2IP (NM_138709) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: DAB2IP
Synonyms: AF9Q34; AIP-1; AIP1; DIP1/2
ACCN: NM_138709
Insert Size: 407 bp
Insert Sequence: >SC201477 3'UTR clone of NM_138709
The sequence shown below is from the reference sequence of NM_138709. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGCGAGTTCAGAAACAGCAGCAATTGTAACTGCCTGAGGAGGGAGGAAGCTACCCAAGGAGAGGGGG
ACTATGGTGGCCAAGGGCAGGGTCTCGGCCTGGGAGGCACCCACGGTTGCAGCCCAGCGGGGTGTC
AGGAGGCCGAGCCTCCCTCCCTGCCGCTGTCCAGGAGCGGCCAGAGGGAGCCACCAGAGACTGAA
GCAGCGTGAGGCGAGGTACCAGCCGCTCCCTGTGGGGTGCGGGCAGAAGAGACTGCACGCTGGGGAGT
GGGGACAGCCTGATGGGGCAGGGGGCCTGCCAAAAATATGTCTGTGTTGTTCTGAATGTGGTGTGTCCT
TGTCTCTGGATCTGGCCGAGTGCATGTGTCCCCCACACCTGTGCCAGGGAGGGGGCTTC
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_138709.2](#)



[View online »](#)

Summary: DAB2IP is a Ras (MIM 190020) GTPase-activating protein (GAP) that acts as a tumor suppressor. The DAB2IP gene is inactivated by methylation in prostate and breast cancers (Yano et al., 2005 [PubMed 15386433]).[supplied by OMIM, May 2010]

Locus ID: 153090

MW: 14.3