

Product datasheet for SC201491

OCIAD2 (NM_152398) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	OCIAD2 (NM_152398) Human 3' UTR Clone
Symbol:	OCIAD2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152398
Insert Size:	542 bp
Insert Sequence:	>SC201491 3'UTR clone of NM_152398 The sequence shown below is from the reference sequence of NM_152398. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
GCCTCCTTACCTGTGAGGAATGCAAAATAAAGCATGGATTAAGTGAGAAGGGAGACTCTCAGCCTTCAG
CTTCCTAAATTCTGTGTCTGTGACTTTTGAAGTTTTTAAACCTCTGAATTTGTACACATTTAAATTT
CAAGTGTACTTTAAATAAAATACTTCTAATGGAACAAAAACATTGTGTCTGTCACTCACTGAAATGTA
TTCCGTGGAAGGCTTTCAAGTGAACAAAGCGTTACATGTTTTCTGCTCTTTAAGAACATGTCGGGGCA
TGGCGCAGTGGCTCATGCCTGTAATCCAGTACTTTGGGAGGCCGAGGCAGGCAGATCATTTGAGGTCA
GGAGTTCATGACCAGCTTGGCCAACATGGTGAAACCCCATCTCTGCTAAAAATATAAAAAATTAGCCAGG
CGTGGTGGCACATGCCTGTAATCCAGCTACTGAGAGGGAGGCTGAGGCAGGAGAATCACTTGAACCCA
GGAGGCGAAAGTTGCAGCGAGCAGAGATTGTGCCACTGCACTCCAGCCTGGGTGCCAGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_152398.4

Locus ID: 132299

MW: 20.5