

## Product datasheet for **SC201500**

### UBL7 (NM\_032907) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	UBL7 (NM_032907) Human 3' UTR Clone
Symbol:	UBL7
Synonyms:	BMSC-UbP; TCBA1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032907
Insert Size:	143 bp
Insert Sequence:	>SC201500 3'UTR clone of NM_032907 The sequence shown below is from the reference sequence of NM_032907. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA <b>GCGATCGCC</b> GAGCTCATCTTTGCTGGAGGAGCCCA <b>TGA</b> ACTCCCTGCTTCCCCTGAACCCCAAGTTGCAGAGG CTACTGCCCTTGGGAGGCACTCATGAAGGTGCCTCCATCTCTCCCTTCCCAATATACCTGATGGTCAA CTCTC <b>ACGCGT</b> AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_032907.5</a></u>
Locus ID:	84993



[View online »](#)

MW: 4.9