

Product datasheet for SC201579

C1orf53 (NM_001024594) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C1orf53 (NM_001024594) Human 3' UTR Clone
Symbol:	C1orf53
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001024594
Insert Size:	172 bp
Insert Sequence:	<p>>SC201579 3'UTR clone of NM_001024594</p> <p>The sequence shown below is from the reference sequence of NM_001024594. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC AAGCAATTCATTCATATTTTATGTTTGACAAGAATTCATCTCTGTTCCCTAACTGTGCTTGTATTT TTAAAAAATAAAGCCCAATTCAGAATTGCTGGATTATTAGTACTTGAACACTAGTTTAATCCTAAAT ACATATATTTAAAGAACATCAATAAAATGAAAAA ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001024594.3</u>
Locus ID:	388722
MW:	6.9


[View online »](#)