

## Product datasheet for SC201591

### SH3TC1 (NM\_018986) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** SH3TC1 (NM\_018986) Human 3' UTR Clone  
**Symbol:** SH3TC1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pMirTarget (PS100062)  
**ACCN:** NM\_018986  
**Insert Size:** 178 bp  
**Insert Sequence:** >SC201591 3'UTR clone of NM\_018986  
 The sequence shown below is from the reference sequence of NM\_018986. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CCCTGGTTGGCCCCCAGCCACCCTCGCGAGGACAGCATCCAAGGGAGTGGGTTTTGTGCAAGGGCTGG
GGGTCTCTGCTCTCTGGTGTGCGCGGTGGCTCATTTTCTGGCAAATGGAGGCACGAACGCAGGGGC
CAAATAGCAATAAATGGGTTTTGTTTTTTTTTGAATAA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

**Restriction Sites:** SgfI-MluI  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).  
**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.  
**RefSeq:** [NM\\_018986.5](#)  
**Locus ID:** 54436  
**MW:** 6.6


[View online »](#)