

Product datasheet for **SC201655**

COPE (NM_199444) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	COPE (NM_199444) Human 3' UTR Clone
Symbol:	COPE
Synonyms:	epsilon-COP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_199444
Insert Size:	195 bp
Insert Sequence:	>SC201655 3'UTR clone of NM_199444 The sequence shown below is from the reference sequence of NM_199444. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC CTGGTGCTACAGTACGCTCCCAGCGCC TGA GGCTGGCCCAGAGCTGTCAGGACCATGAAGCCAGGACAG AGGCCAGGAGCCAGCCCTGCAGCCCTCCCCACCCGGCATCCACCTGCATCCCCTCTGGGGGAGGAGCC CACCCCAGCACCCCATCTGTTAATAAATATCTCAACTCCAGGGTGTTCACCTGA ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_199444.2</u>



[View online »](#)

Summary: The product of this gene is an epsilon subunit of coatomer protein complex. Coatomer is a cytosolic protein complex that binds to dilysine motifs and reversibly associates with Golgi non-clathrin-coated vesicles. It is required for budding from Golgi membranes, and is essential for the retrograde Golgi-to-ER transport of dilysine-tagged proteins. Coatomer complex consists of at least the alpha, beta, beta', gamma, delta, epsilon and zeta subunits. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Locus ID: 11316

MW: 6.6