

Product datasheet for SC201696

C2orf81 (NM_001145054) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C2orf81 (NM_001145054) Human 3' UTR Clone
Symbol:	C2orf81
Synonyms:	hCG40743
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001145054
Insert Size:	182 bp
Insert Sequence:	<p>>SC201696 3'UTR clone of NM_001145054</p> <p>The sequence shown below is from the reference sequence of NM_001145054. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC ATCCAGACAGGTGCCCCAAAGCCCAGGTACCGTAGCACAGCTAATGAAGAACTCAGCCCCCAAAGTGT GGTCACGCTCCTCTAAGCCTGCTGCCTCCCTCTGATCCCTGAGCCTCTGGCAGTAGTAAGTGGTCCCTC CCTCTGCTAGCCAGAAATAAACACCTGAGTTGCCTTAGGAACTA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001145054.2</u>
Locus ID:	388963


[View online »](#)

MW: 6.5