

Product datasheet for **SC201711**

HPR (NM_020995) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	HPR (NM_020995) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	HPR
Synonyms:	A-259H10.2; HP
ACCN:	NM_020995
Insert Size:	195 bp
Insert Sequence:	>SC201711 3'UTR clone of NM_020995 The sequence shown below is from the reference sequence of NM_020995. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TGGTTTCAGAAGACCATAGCTGAGAACTAATGCAAGGCTGGCCGGAAGCCCTTGCCCTGAAAGCAAGATT TCAGCCTGGAAGAGGGCAAAGTGGACGGGAGTGGACAGGAGTGGATGCGATAAGATGTGGTTTGAAGCT GATGGGTGCCAGCCCTGCATTGCTGAGTCAATCAATAAAGAGCTTCTTTTGACCCA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_020995.4



[View online »](#)

Summary: This gene encodes a haptoglobin-related protein that binds hemoglobin as efficiently as haptoglobin. Unlike haptoglobin, plasma concentration of this protein is unaffected in patients with sickle cell anemia and extensive intravascular hemolysis, suggesting a difference in binding between haptoglobin-hemoglobin and haptoglobin-related protein-hemoglobin complexes to CD163, the hemoglobin scavenger receptor. This protein may also be a clinically important predictor of recurrence of breast cancer. [provided by RefSeq, Oct 2011]

Locus ID: 3250

MW: 7.4