

Product datasheet for SC201719

CCHCR1 (NM_001105563) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCHCR1 (NM_001105563) Human 3' UTR Clone
Symbol:	CCHCR1
Synonyms:	C6orf18; HCR; pg8; SBP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001105563
Insert Size:	183 bp
Insert Sequence:	<p>>SC201719 3'UTR clone of NM_001105563</p> <p>The sequence shown below is from the reference sequence of NM_001105563. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC TCCAGCTCCAATCCCAGATGAGCAGCTAAGCAGCTGACAGTTGGAGGGAAAGCCAGCCTGGGGGCTGG GAGGATCCTGGAGAAGTGGGTGGGGACAGACCAGCCCTTCCCATCCTGGGGTTGCCCTGGGGGATACC AGCTGAGTCTGAATTCTGCTCTAAATAAGACGACTACAGAAGGA ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001105563.3


[View online »](#)

Summary:	This gene encodes a protein with five coiled-coil alpha-helical rod domains that is thought to act as a regulator of mRNA metabolism through its interaction with mRNA-decapping protein 4. It localizes to P-bodies, the site of mRNA metabolism, with an N-terminus that is required for this subcellular localization, suggesting it is a P-body component. Naturally occurring mutations in this gene are associated with psoriasis. [provided by RefSeq, May 2017]
Locus ID:	54535
MW:	6.3