

Product datasheet for **SC201742**

IFI35 (NM_005533) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: IFI35

Synonyms: IFP35

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_005533

Insert Size: 187 bp

Insert Sequence: >SC201742 3'UTR clone of NM_005533
The sequence shown below is from the reference sequence of NM_005533. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTAGCAGTCTTCACCTCTGAGTCAGGCTAGGGGCTCCCTTCTCATCCTCCCCACCCCCCGCCAAGG
TTCTCACACTGGCCTGGGCTTGGGTGCCCATATAGGAGGTCTGTATGTTACCAACAGTGCGGAGGGGT
CACACATTGCAAAACACTGCCGAGAACAGTAAAAAGAGCCTGCATGCCA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_005533.5
Summary:	Not yet known.[UniProtKB/Swiss-Prot Function]
Locus ID:	3430
MW:	6.7