

Product datasheet for **SC201761**

SMIM29 (NM_001008703) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SMIM29 (NM_001008703) Human 3' UTR Clone
Symbol:	SMIM29
Synonyms:	C6orf1; LBH
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001008703
Insert Size:	438 bp
Insert Sequence:	>SC201761 3'UTR clone of NM_001008703 The sequence shown below is from the reference sequence of NM_001008703. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC CGGATGCCACTGCTGGATGTCAAGACGTGACCTGACCCCTTGCCCCACCCTTCAGAGCCTGGGGTCT GGACTGCCTGGGGCCCTGCCATCTGCTTCCCCTGCTGTACCTGGCTCCCCCTGCTGGGTGCTGGGTCT CCATTTCTCCCTCACCCACCCTCAGCAGCATCTGCTTCCCATGCCCTCACCATCACCTCACTGCCCCC AGGCCTTCTGCCCTTTGTGGGTGTTGAGCTCACCGCCACCCACAGGCACTCATGGGAAGAGGCTTTCC TTCTGGGATGGGGCGGCTGGTAGACACCTTTGCTTTCTCTAGCCCTCCTGGGCTGGGCTTGGGCACAA ATCCCCAGGCAGGCTTTGGAGTTGTTTCCATGGTGATGGGGCCAGATGTATAGTATTCAGTATATATTT TGTAATAAAATGTTTTGTGGCTA ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_001008703.4

Locus ID: 221491

MW: 16