

Product datasheet for **SC201802**

MTH1 (NUDT1) (NM_198954) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MTH1 (NUDT1) (NM_198954) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	NUDT1
Synonyms:	MTH1
ACCN:	NM_198954
Insert Size:	204 bp
Insert Sequence:	>SC201802 3'UTR clone of NM_198954 The sequence shown below is from the reference sequence of NM_198954. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TACACACTCCGCGAGGTGGACACGGTCTAGCGGGAGCCCAGGGCAGCCCCTGGGCAGGAGACGTGGCTG CTGAACAGCCGCAAACCATCTTACCTGGGGGCATTGAGTGGCGCAGAGCCGGTTTCATCTGGAATTA ACTGGATGGAAGGAAAAATAAAGCTATCTAGCGGTGAAAAAAAAAAAAAAAAAAAAAAAAAAAAA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_198954.1</u>



[View online »](#)

Summary:

Misincorporation of oxidized nucleoside triphosphates into DNA/RNA during replication and transcription can cause mutations that may result in carcinogenesis or neurodegeneration. The protein encoded by this gene is an enzyme that hydrolyzes oxidized purine nucleoside triphosphates, such as 8-oxo-dGTP, 8-oxo-dATP, 2-hydroxy-dATP, and 2-hydroxy rATP, to monophosphates, thereby preventing misincorporation. The encoded protein is localized mainly in the cytoplasm, with some in the mitochondria, suggesting that it is involved in the sanitization of nucleotide pools both for nuclear and mitochondrial genomes. Several alternatively spliced transcript variants, some of which encode distinct isoforms, have been identified. Additional variants have been observed, but their full-length natures have not been determined. A rare single-nucleotide polymorphism that results in the production of an additional, longer isoform (p26) has been described. [provided by RefSeq, Dec 2018]

Locus ID:

4521

MW:

7.9