

Product datasheet for **SC201819**

DOCK9 (NM_001130049) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DOCK9 (NM_001130049) Human 3' UTR Clone
Symbol:	DOCK9
Synonyms:	ZIZ1; ZIZIMIN1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001130049
Insert Size:	189 bp
Insert Sequence:	>SC201819 3'UTR clone of NM_001130049 The sequence shown below is from the reference sequence of NM_001130049. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC ATCACTCTGAAAGTGTCTTGGAGCCAA TAG TTGGTGAACGTGTCACACTTGTGTGGTAGGACCTTGAAG TCTAAGTTGCTTTCCTGAGTATTCTTTTCCTGCTTGTGATAGTCAACAACCTGAAACCCCTCAGCCATGC CCTGAAATAAAGGTCCCGGATGCCTGTGACTCCTCAGGATCATAAGTTAA ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001130049.2</u>



[View online »](#)

Summary: Guanine nucleotide-exchange factor (GEF) that activates CDC42 by exchanging bound GDP for free GTP. Overexpression induces filopodia formation.[UniProtKB/Swiss-Prot Function]

Locus ID: 23348

MW: 6.8