

Product datasheet for **SC201860**

REX1BD (NM_001100419) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	REX1BD (NM_001100419) Human 3' UTR Clone
Symbol:	REX1BD
Synonyms:	C19orf60
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001100419
Insert Size:	187 bp
Insert Sequence:	>SC201860 3'UTR clone of NM_001100419 The sequence shown below is from the reference sequence of NM_001100419. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC CTTAGGTTTGTATGCGGAATCTGCCGAG TGAT GGCGGCTCCCCAGGGATGCGCCGAGGGAGATGGGAAAC GGGGCGGATGGCGCCAGCCAGCCCTAACTGCCAGCTGGCTGGGGTTGCGCCCCACTGCGCTGCTGAC CTTCTGCAGTTCAGACACCTCCACAATAAAGAGCTCCTCCTCTGTA ACGCGT AAGCGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001100419.2</u>
Locus ID:	55049



[View online >](#)

MW: 6.6