

Product datasheet for **SC201868**

SCRN2 (NM_138355) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SCRN2 (NM_138355) Human 3' UTR Clone
Symbol:	SCRN2
Synonyms:	Ses2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_138355
Insert Size:	182 bp
Insert Sequence:	>SC201868 3'UTR clone of NM_138355 The sequence shown below is from the reference sequence of NM_138355. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC GTGAAGAGGGAGAGCCAGGCTTATGCG TA AGCTTCATAGCTTCTGCTGGCCTGGGGTGGACCCAGGACC CCTGGGGCCTGGGTGCCCTGAGTGGTAAAGTGGAGCAATCCCTTACGCTCCTTGCCATGTTCTG AGCGGCCAGCTTGGCCTTTGCCTAATAAATGTGCTTTATTTTC ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_138355.4</u>
Locus ID:	90507



[View online »](#)

MW: 6.5