

Product datasheet for SC202237

HEATR7A (MROH1) (NM_032450) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: HEATR7A (MROH1) (NM_032450) Human 3' UTR Clone
Symbol: HEATR7A
Synonyms: HEATR7A
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_032450
Insert Size: 216 bp
Insert Sequence: >SC202237 3'UTR clone of NM_032450
 The sequence shown below is from the reference sequence of NM_032450. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GCCCTGGGCCGCTGGTGAAGCTCGCTTAGGCTCCGGCCAGCACCCCCAGCGAGGAGTATGCACCCAG
ACCTGTGCCTGAGCTCCAAGACAGGGCCTCCTGAGGACCACAGCCTGGGCACACGACTGGAGGGGCCTG
GCCCCAGAACAGGCACTGCTGGGGACCAACCCAAGCCCTTCAGTGAGGGATGTGCCCAATAAACTCAT
CTGCTCCCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq: NM_032450.3
Locus ID: 727957


[View online »](#)

MW: 7.7