

Product datasheet for **SC202299**

REXO5 (NM_001144924) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: REXO5

Synonyms: NEF-sp

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_001144924

Insert Size: 215 bp

Insert Sequence: >SC202299 3'UTR clone of NM_001144924
The sequence shown below is from the reference sequence of NM_001144924. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAAAGCGCTGGCCAGGCCTGTGTTCTGAGTCGGCCTGCCATGTTTCCATGTGCCATTTCTTACCCCT
TGTAAGCAATGGCAAAGAATGTGGTCAGGCTGTAGCCTCCCAACCAGCAGACAGCTTTATGGAACTT
GGTATAGCAGCTAAAAGAGTTTAGTTTGTATATGGCATGTATAAGTTTCAATAAATGCCTAAAGTT
CAAGCATA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001144924.2](#)

Locus ID: 81691

MW: 7.9