

Product datasheet for **SC202301**

REXO5 (NM_030941) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	REXO5
Synonyms:	NEF-sp
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_030941
Insert Size:	215 bp
Insert Sequence:	>SC202301 3'UTR clone of NM_030941 The sequence shown below is from the reference sequence of NM_030941. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GAAAGCGCTGGCCAGGCCTGTGTTCTGAGTCGGCCTGCCATGTTTCCATGTGCCATTTCTTACCCCT TGTAGGCAATGGCAAAGAATGTGGTCAGGCTGTAGCCTCCCAACCAGCAGACAGCTTTATGGAACTT GGTATAGCAGCTAAAAGAGTTTAGTTTGTATATGGCATGTATAAGTTTCAATAAATGCCTAAAGTT CAAGCATA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_030941.3](#)

Locus ID: 81691

MW: 7.9