

Product datasheet for **SC202476**

APOA1BP (NAXE) (NM_144772) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	APOA1BP (NAXE) (NM_144772) Human 3' UTR Clone
Symbol:	APOA1BP
Synonyms:	AIBP; APOA1BP; PEBEL; YJEFN1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_144772
Insert Size:	245 bp
Insert Sequence:	>SC202476 3'UTR clone of NM_144772 The sequence shown below is from the reference sequence of NM_144772. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC GACACCGAGTGTGTCTATCGTCTGCAG TGA GGGAAGGTGGGTGGTATTCTTCCCAATAAAGACTTAGA GCCCTCTCTCCAGAACTGTGGATTCTGGGAGCTCCTCTGGCAATAAAAGTCAGTGAATGGTGGAAAG TCAGAGACCAACCCTGGGATTGGGTGCCATCTCTCTAGGGGTAAACACAAAGGGCAAGAGGTTGCTATG GTATTTGGAAACAATGAAAATGGACTGTTAGATGCCAA ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_144772.3</u>



[View online »](#)

Summary: The product of this gene interacts with apolipoprotein A-I (apoA-I), the major apolipoprotein of high-density lipoproteins (HDLs). It is secreted into some bodily fluids, and its synthesis and secretion are stimulated in vitro by incubating cells with apoA-I. The human genome contains related pseudogenes. [provided by RefSeq, Jul 2008]

Locus ID: 128240

MW: 8.9