

Product datasheet for **SC202576**

POLR1H (NM_014596) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	POLR1H (NM_014596) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	POLR1H
Synonyms:	A12.2; HTEX-6; HTEX6; hZR14; Rpa12; tctex-6; TCTEX6; TEX6; ZNRD1; ZR14
ACCN:	NM_014596
Insert Size:	240 bp
Insert Sequence:	>SC202576 3'UTR clone of NM_014596 The sequence shown below is from the reference sequence of NM_014596. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TGCAAGTTCCAGGAGAAGGAAGACTCT TGAC CTTTTTCTGGCAACTCTACAGTCCCTCCCTCCTTTC GGAAGGTGAAGGATACTGGGTTTTAGATGCCTTGCCATCCTGTCTGGTTGCAATGTTTTGCTCCCAG AAGAGAATCAGATCATCATGTGGGATTACCATTGTTCTGGAGTACTCCTACCCTTAGTTGAATTTCC TTATTAAGTTATATTTTTCTATAAGACCCTGA ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_014596.6</u>



[View online »](#)

Summary: This gene encodes a DNA-directed RNA polymerase I subunit. The encoded protein contains two potential zinc-binding motifs and may play a role in regulation of cell proliferation. The encoded protein may be involved in cancer and human immunodeficiency virus progression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]

Locus ID: 30834

MW: 9.1