

## Product datasheet for **SC202705**

### **KHDC1 (NM\_030568) Human 3' UTR Clone**

#### Product data:

Product Type:	3' UTR Clones
Product Name:	KHDC1 (NM_030568) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	KHDC1
Synonyms:	bA257K9.4; C6orf147; C6orf148; Em:AC019205.8; KHDC1L; NDG1
ACCN:	NM_030568
Insert Size:	245 bp
Insert Sequence:	>SC202705 3'UTR clone of NM_030568 The sequence shown below is from the reference sequence of NM_030568. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GGTTCCATTTGAGTTCCTGTATCCGTAATTAAGGACTGATCTCCCTGGGAAGTGGGGATGAGAGGA GAAATCTTGGCTTCACTTTTGGTGAGAAAAGACATGATAGACACTGAGGAAAGTTGAACTTGTGGCCGG GACCTGTCCCTCTCAGGTCTTCTAGACATGTCTGAGGGAAGTCTTGAATGTAAGAAAGTGTATTGG CAAAGCAATCCAATAAAAATAAGCTACCTATGATGCAAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_030568.5</a></u>
Locus ID:	80759
MW:	9.1



[View online »](#)