

Product datasheet for **SC202763**

CRELD2 (NM_024324) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CRELD2 (NM_024324) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CRELD2
ACCN:	NM_024324
Insert Size:	260 bp
Insert Sequence:	>SC202763 3'UTR clone of NM_024324 The sequence shown below is from the reference sequence of NM_024324. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_024324.5
Summary:	May regulate transport of alpha4-beta2 neuronal acetylcholine receptor.[UniProtKB/Swiss-Prot Function]
Locus ID:	79174
MW:	9.6



[View online »](#)