

Product datasheet for **SC202968**

CTAGE4 (NM_198495) Human 3' UTR Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | 3' UTR Clones |
| Product Name: | CTAGE4 (NM_198495) Human 3' UTR Clone |
| Symbol: | CTAGE4 |
| Synonyms: | cTAGE-4 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pMirTarget (PS100062) |
| ACCN: | NM_198495 |
| Insert Size: | 273 bp |
| Insert Sequence: | >SC202968 3'UTR clone of NM_198495 The sequence shown below is from the reference sequence of NM_198495. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site |

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CCTGGATTTTACCCCAACCCACATTCTGAAGGTAGAAGCGAGTTCCTTCAGGATTGATCCGCCTTC
AAAGGAGCCTGCTACTGGACATCCAGAACCACAGCAAGAAACCTGACAATATTGTTGCTTTCTTCAAAA
GTAATTTTACTGATCTCATTTTCAGTTAAGTAACTGCTGTTACTTAAGTGATTGCACTTTTCTCAAA
TTGAAGTTTAATGGAATAATAGTTCTCAGGATAGTATTTTCGTAATAAAGATGGTTTGAATATGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

| | |
|--------------------|---|
| Restriction Sites: | Sgfl-MluI |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | NM_198495.3 |
| Locus ID: | 100128553 |



MW: 10.7