

## Product datasheet for **SC203021**

### CD64 (FCGR1B) (NM\_001004340) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	CD64 (FCGR1B) (NM_001004340) Human 3' UTR Clone
Symbol:	CD64
Synonyms:	CD64b; FCG1; Fcgamma1a; FCGR1; FCGR1A; FcRI; IGFR1; IGFRB
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001004340
Insert Size:	1039 bp
Insert Sequence:	>SC203021 3'UTR clone of NM_001004340 The sequence shown below is from the reference sequence of NM_001004340. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CACCGGAAGGAGCCCCAGGGGGCCACGTAGCAGCGGCTCAGTTGGTGGCCATCGATCTGGACCGTCCCC
TGCCCACTTGCTCCCGTGAGCACTGCGTACAAACATCCAAAAGTTCAACAACACCAGAAGTGTGTGC
TCATGGTATATACTCTTAAAGCAAATAAATGAAGTACTTCAACTGGGATACATTTGAAATGTGGTC
ATCAAAGATGACTTGAAATGAGGCTACTCTAAAGAATTCTTAAAAACTTACAAGTCAAGCCTAGCCT
GATAATCCTATTACATAGTTTAAAAATAGTATTTTATTTCTCAGAACAAAGGTAAAAAGGTGAGTGGGT
GCATATGTACAGAAGATTAAGACAGAGAAACAGACAGAAAGAGACACACACAGCCAGGAGTGGGTAG
ATTTACAGGAGACAAGAGGGAATAGTATAGACAATAAGGAAGGAAATAGTACTTACAAATGACTCCTAA
GGGACTGTGAGACTGAGAGGGCTCACGCCTCTGTGTTCCAGGATACTTAGTTTCATGGCTTTTCTTTGA
CTTTACTAAAAGAGAATGTCTCCATACGCGTCTAGGCATACAAGGGGTAAGTCAATGATGAGAAATGG
ATGTGTTATTCTTGCCCTCTCTTTGAGGCTCTCTATAACCCCTCTATTCTAGAGACAACAAAAATG
TTGCCAGTCTCAGGCCCTGCCCTGTAGGAAGGCAGAATGTAAGTCTTTTGTAAACGATTAAGT
CCAAATCTCCAAGTCCGCACTGCAAAGAGACGCTTCAAGTGGGAGAAGCGCGATATCATAGAGTCC
AGATCTTGCCCTCAGAGATTTGCTTTACCTTCTGATTTTCTGGTACTAATTAGCTTCAGGATACGCT
GCTCTCATACTTGGGCTGTAGTTTGGAGACAAAATATTTTCTGCCACTGTGTAACATAGCTGAGGTAA
AACTGAACTATGTAATGACTCTACTAAAAGTTTAGGGAAAAAAACAGGAGGAGTATGACACACACA
GCAA
AGCGGACCGACTTACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCC
CAACCTGCCATCACGAGATTTGATTCCACCGCCGC

```

Restriction Sites: SgfI-RsrII



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_001004340.4</a>
<b>Summary:</b>	Three distinct, but closely related classes of receptors that bind the Fc portion of IgG have been identified (Fcgamma RI, II and III). The FcgammaRI family consists of three closely related genes termed A, B, and C. This gene likely encodes a non-functional protein that is not detectable at the cell surface and binds ligand with low affinity. [provided by RefSeq, Nov 2019]
<b>Locus ID:</b>	2210
<b>MW:</b>	39.2